


Outside the cytoplasm: Biophysical ligand screening for membrane proteins?

September 2010
Gregg Siegal

Leiden University & ZoBio
Leiden, The Netherlands




Universiteit Leiden


Gregg Siegal
ZoBio



Fragments: Why small is beautiful




Why biophysical fragment screening for GPCRs?




Gregg Siegal
ZoBio

The TINS method for finding hits


TINS = Target Immobilized NMR Screening



The TINS Screening Station




Dual cell
Sample holder



Gregg Siegal
ZoBio

An example of raw data



Summary of Selected Immobilized Targets


Protein	Size	Functional	TINS hits
Protease	44 kDa	LB	5.9%
HSP90	24 kDa	LB	6.5%
Small GTPase	20 kDa	LB/BA	9% apo form 3% NDP form
Viral enzyme	67 kDa	LB/BA	9.5%
DsbB (Bacterial mem. Prot.)	14 kDa	BA	7.3%
Various kinases (pY,pS/T)	30-35 kDa	LB/BA	3.8-5.1%
KcsA (Ion channel)	57 kDa	LB	Feasibility only (95 cmpds, 7%)
Metalloproteins	105 kDa homotrimer	LB/BA	5-8.5%
Prot-Prot Interaction (5)	14-100 kDa	LB	3-6%

LB ligand binding
BA biochemical assay


Gregg Siegal
ZoBio



Fragment screening with TINS



Undruggable Target
11 hits





Druggable Target
89 hits

Methods for Hit Validation

- Affinity ranking


- Competition Binding
- SPR
- HSQC Binding Site Determination

Hit Validation: Competition Binding in TINS





Hit Validation: SPR (Biacore T200)

Biacore Analysis of two small molecules binding to MDM4





Hit Validation: K_D and Binding Site Characterization by HSQC

- Chemical shift perturbation for each HN
 - magnitude of resonance shift due to compound binding
- $CSP = \text{sqrt} [\Delta H_{\text{ppm}}^2 + (\Delta N_{\text{ppm}}/6.5)^2]$



CSP observation on 80 TINS Hits

Low resolution binding site determination



CSP: ■ > ■ > ■

Hits bind in the vicinity of PPI site



Gregg Siegal
ZoBio



The Challenges of Fragment Discovery on MPs

- Protein production/solubilization/stability
- Non-specific binding
- Slow kinetics
- Fragment size

FBDD on Membrane Proteins: DsbB



Immobilized DsbB has 90% the activity of the soluble protein.

Inaba et al., Cell, 2006, **127**, p.789

Target: E. coli DsbB

Ref: E.coli OmpA

Both proteins DPC solubilized

Cmpds Screened: 1,270

Protein used: 2 mg



Hits: 93

Früh et al, Chem. Biol., 2010, v. 17, p. 881

Gregg Siegal
ZoBio




TINS Ligand Screening on Membrane Proteins




1,270 fragments tested → 7.3 % hit rate → Validation/characterization


Membrane Protein Inhibitors: Mode of Action




Frag	IC50 (uM)	Hill Slope
1	7 ± 1	0.80 ± 0.10
2	10 ± 1	0.80 ± 0.10
3	50 ± 10	0.80 ± 0.10
4	70 ± 10	1.00 ± 0.10
5	100 ± 10	1.40 ± 0.10
6	115 ± 15	1.15 ± 0.05
7	170 ± 10	1.40 ± 0.10
8	190 ± 10	1.20 ± 0.10







Fragment Number	Fragment [Conc]	Substrate	K_{cat} (M DsbA/ M DsbB-s ⁻¹)	K_m (μM)
2	0 μM	UQ1	3.7 ± 0.1	1.6 ± 0.1
	75 μM	UQ1	3.1 ± 0.2	13.2 ± 1.7
2	0 μM	DsbA	4.3 ± 0.1	2.4 ± 0.2
	75 μM	DsbA	2.5 ± 0.1	2.8 ± 0.1
8	0 μM	UQ1	4.0 ± 0.1	2.2 ± 0.0
	75 μM	UQ1	2.3 ± 0.1	9.7 ± 1.5
8	0 μM	DsbA	4.3 ± 0.1	1.9 ± 0.1
	75 μM	DsbA	2.4 ± 0.2	4.1 ± 0.5







Structural Model of Fragment Binding






Competitive






Competitive




Mixed Model




Effects of the micelle on screening



Gregg Siegal
ZoBio



Effects of the micelle on screening




	Obs/Unobs (cLogP)	Hits	cLogP	BioAssay in ND	BioAssay in DPC
Micelle	127/56 (0.9/1.8)	8	1.34	-	+
NanoDisc	164/19 (1.1/1.6)	8	2.21	++	++
Both	-	14	2.13	++	++

Gregg Siegal
ZoBio




Fragment Screening of GPCRs: Adenosine 2a & β 1 Adrenergic receptors



$A_{2a}R$ agonists – anti-inflammatory therapeutic potential

$A_{2a}R$ antagonists – used to treat Parkinson's disease as $A_{2a}R$ dimerises with dopamine D_2 receptor

Radioligand binding on immobilized adenosine A_{2a} receptor




Immobilized hA2aR maintained 60% the activity in five days.


Profile of Ligand Binding in the Screen

T/R distribution by compound


531 fragments assayed against A2aR
94 hits



Gregg Siegal
ZoBio




Hit validation by equilibrium radioligand displacement



Each hit assayed at 500 μ M for 3 H- ZM241385 displacement using wild type A2aR.


ID	K_i (M)	EC_{50} (M)	Hill Slope	T/R ratio
ZB643	$3.0 \pm 0.1 \times 10^{-6}$	$7.0 \pm 0.3 \times 10^{-6}$	-0.5 ± 0.1	0.65
ZB418	$3.0 \pm 0.1 \times 10^{-5}$	$5.5 \pm 0.1 \times 10^{-4}$	-2.0 ± 0.5	0.44
ZB1703	$4.1 \pm 0.1 \times 10^{-5}$	$6.4 \pm 0.1 \times 10^{-4}$	-1.1 ± 0.4	0.60
ZB1166	$1.2 \pm 0.1 \times 10^{-4}$	$1.5 \pm 0.1 \times 10^{-3}$	-0.9 ± 0.1	0.61
ZB2086	$1.2 \pm 0.2 \times 10^{-4}$	$2.4 \pm 0.3 \times 10^{-3}$	-0.9 ± 0.3	0.74
ZB114	$8.2 \pm 0.3 \times 10^{-5}$	$1.2 \pm 0.1 \times 10^{-3}$	-2.2 ± 5.4	0.59
ZB1605	$8.7 \pm 0.3 \times 10^{-5}$	$1.4 \pm 0.3 \times 10^{-3}$	-1.1 ± 0.5	0.37
ZB1967	$3.2 \pm 0.2 \times 10^{-4}$	$6.1 \pm 0.1 \times 10^{-3}$	-0.8 ± 1.3	0.61

Mode of action: Orthosteric vs Allosteric modulators




Screen all hits at $t = 50\%$ ZM bound.

TINS hits as negative allosteric modulators (NAMs) of A2aR




	ZM241385	ZM241385 + 2.5mM ZB2044	ZM241385 + 2.5mM ZB418
k_{off} of ZM241385 (min^{-1})	0.010 ± 0.003	0.025 ± 0.004	0.048 ± 0.011
Half Life (min)	69	27	15

TINS hits as positive allosteric modulators (PAMs) of A2aR




	ZM241385	ZM241385 + 2.5mM ZB178	ZM241385 + 2.5mM ZB728
k_{off} of ZM241385 (min^{-1})	0.0100 ± 0.0026	0.0034 ± 0.0016	0.0057 ± 0.0023
Half Life (min)	69	205	123

Structural information from TINS?



$\beta 1$ AR



Gregg Siegal
ZoBio

Summary

- True fragment screening (e.g. biophysical) can be successfully applied to GPCRs.
- Novel orthosteric ligands can be found using FBDD.
- Novel allosteric modulators of GPCR function can be found using FBDD.

Applications

- Primary screen for novel matter
- Scaffold hopping
- Allosteric modulators

Acknowledgements

DsbB

Leiden University

Virginie Früh

University of Virginia

Yunpeng Zhou

John Bushweller

ZoBio

Eiso AB

β 1AR/A2aR

Leiden University

Francis Figaroa Laura Heitman, Ad Ijzerman

ZoBio

Dan Chen, Johan Hollander, Eiso AB

MRC Laboratory of Molecular Biology

Chris Tate, Tony Warne, Maria Serano-Vega

Heptares Therapeutics

Markus Koglin, James Errey, Andrei Zhukov, Fiona Marshall

NanoDiscs

University of Illinois-Urbana Champagne

Yelena Grinkova
Stephen Sligar



Gregg Siegal
ZoBio

